


Mini Extruder Pro

Extrusion technique guide



Avanti Research
A CRODA BRAND

A quick-reference workflow for preparing liposomal suspensions with confidence.



Membrane selection matters

Extrusion of multilamellar liposomal suspensions using membranes with a pore size $>0.2\mu\text{m}$ does not produce unilamellar liposomes. Liposomes produced with larger pore membranes will yield a polydispersed suspension of multilamellar liposomes. Unilamellar liposomal suspensions with a low polydispersity can only be prepared with membranes having a pore size of $\leq 0.2\mu\text{m}$.

Pre-wet first

Pass a syringe full of buffer through the extruder, then discard the buffer

Single use

Polycarbonate membranes and filter supports are for one liposome preparation only.

80°C max

Do not exceed 80°C; higher temperatures can damage the extruder and syringe components.

Setup workflow

1 Prepare dry lipid mixture

Prepare by lyophilization or evaporation

2 Heat the base

Place base on hot plate and allow block to reach temperature (~15 min).

3 Hydrate sample

Hydrate at least 30 mins above phase transition temperature.

4 Load and dock syringes

Once hydrated, load sample into one gas-tight syringe and dock into one side of the Mini Extruder Pro.

Set the empty syringe plunger to zero and dock it on the opposite side.



Do not twist syringes



Slide into body until firmly seated. Twisting can loosen needles and create leaks.

Run the extrusion

Use gentle, controlled transfers and keep the system seated correctly throughout the workflow.

1. Setup + thermal contact

Seat the extruder apparatus into the base and then insert the syringes.
Use the swing arm clips to keep the syringes aligned and in good thermal contact.




2. Seat syringes correctly

Do not twist; syringes should fit tight into the extruder.



3. Equilibrate temperature

Allow suspension to equilibrate for approximately 5-10 minutes.



4. Transfer back and forth

Gently push the filled syringe to transfer the lipid solution to the alternate syringe, then transfer it back.

Minimum 11 passes

More passes improve homogeneity.



5. Collect from the alternate syringe

Collect final extrusion from the alternate syringe. Remove the unit from the heating block, then pull syringes straight out by sliding on the base. Use care: solution and syringe may be hot.

Scan here for full technique guide:



Non-warranty

The information in this publication is believed to be accurate and is given in good faith, but no representation or warranty as to its completeness or accuracy is made. Suggestions for uses or applications are only opinions. Users are responsible for determining the suitability of these products for their own particular purpose. No representation or warranty, expressed or implied, is made with respect to information or products including, without limitation, warranties of merchantability, fitness for a particular purpose, non-infringement of any third-party patent or other intellectual property rights including, without limit, copyright, trademark and designs. Unless otherwise stated, any trademarks identified herein are trademarks of the Croda group of companies.
©2026 Croda International.