


Now Available from Avanti: Amphiphilic Polymethacrylate (PMA) Copolymers



Advantages of PMA Based Lipid Nanodiscs:

Spontaneously Forms Lipid Nanodiscs through Fragmentation of Lipid Bilayer

Amenable to Biophysical Characterization using CD, UV/vis, and Fluorescence Spectroscopy

Ability to Produce Nanodiscs of Various Sizes

Larger PMA Polymer Nanodiscs can be aligned in a magnetic field allowing for solid-state NMR studies on membrane proteins or peptides.

